


## MEMO: How to calculate the limit of detection (LOD)?

- a) **Definition:** the limit of Detection  $LOD$  with a power  $(1 - \beta)$  is defined (in concentration unit) as the minimum concentration for which detecting the target gene in a well is possible with a probability of  $1 - \beta$ .

In other words, this is the minimum concentration that can be said to be non-zero and statistically higher than the limit of blank  $LOB$  with a  $1 - \beta$  probability (typically 95% for  $\beta = 5\%$ ).



**Definition of the sampling limit  $LOS$**  with a power of  $(1 - \beta) = 95\%$ : the lowest concentration in the well for which the probability of having no target gene in the analyzed sub-sample is  $\beta = 5\%$ .

### b) A calculation method:

- 1) Determine the limit of blank  $LOB(95\%)$  of your experiment by following the method of the memo "How to calculate the limit of blank".

- 2) Calculate 
$$p_0 = \frac{2b+z^2+z\sqrt{z^2+4b(1-b/N)}}{2N(1+z^2/N)}$$

where:

- $b = LOB(95\%)$  is the 95% limit of blank
- $z = 1.645$  is the "one-tail" quantile at 95%
- $N$  is the total number of partitions that are generated on average in a well
- $p_0$  is the higher-value solution of the following equation (which can be simplified as a second degree equation in  $p$ ):

$$p = b/N + z \sqrt{p(1-p)/N}$$

3) The LOD with a 95% confidence level is determined as follows (calculation based on the Normal Law approximation and the Poisson Law):

- In concentration in the well (cp / uL):

$$LOD(95\%) = -N \ln(1 - p_0)/V$$

$$LOD(95\%) = -\ln(1 - p_0)/v$$

where  $V = N v$  is the analyzed volume (cumulative total volume of all the partitions generated in the well, with  $v$  the volume of each partition)

- In number of copies included in the volume analyzed in the well:

$$LOD_{cp}(95\%) = [-N \ln(1 - p_0)]$$

knowing that this unit in number of copies can be approximated as "number of positive partitions" in the well.

- For example: if  $LOB(95\%) = 2$ ,  $v = 0.000442$  uL and  $N = 28000$ , then  $LOD(95\%) = 0.49$  cp / uL and  $LOD_{cp}(95\%) = 7$  copies (or "positive partitions")
- Special case: if  $LOB(95\%) = 0$  (i.e. false positives are never present) then  $LOD(95\%) = 3/V$  and  $LOD_{cp}(95\%) = 3$  copies. In fact, if there are never false positives, then the 95% detection limit is equal to the 95% sampling limit which is:  $LOS(95\%) = 3/V$  cp / uL.